

# MARLEN UltraLite™

## One-Piece Disposable, Drainable System For Ileostomies and Colostomies

*Choose the Style, Size,  
Stoma Opening, and Skin Barrier  
That Suits Your Requirement*



Now Available with Marlen's  
Kwick-Klose™ Fastener

*Slimmer Design  
Unparalleled Flexibility  
The Ultimate  
in Comfort*

*Available with either:*

**MARLEN'S AquaTack™ Hydrocolloid Barrier**

- Absorbs skin moisture
- Adheres to moist skin
- Seals securely; removes easily from sensitive skin (controlled adhesion)


**MARLEN'S Skin Shield™ Adhesive Skin Barrier**

- Stays soft
- Elastic and comfortable
- Will not absorb drainages
- Offers long-lasting protection

*or*

*Choose Deep or  
Shallow Convexity*

*Choose Flat, Pre-Cut or  
Cut-To-Fit*




- Cloth-like comfort cover on body side
- Choice of transparent or opaque outer side
- Stoma sizes from 1/2 inch (12mm) to 1-3/4 inches (44mm) or Cut-To-Fit (flat only)
- Completely odor-proof material

# For Stoma Openings Up to 3" **MARLEN Ultra™**

## MARLEN UltraExtra™ One-Piece Disposable System For Ileostomies and Colostomies


- With Marlen's Skin Shield™ Adhesive Skin Barrier
- Stoma sizes from 1-7/8" (48mm) to 3" (76mm) and oval sizes from 1" x 1-1/2" (25mm x 38mm) to 2" x 3" (50mm x 76mm) OR Cut-To-Fit
- Choose small or medium pouch size
- Choice of transparent or opaque outer side



## For Stoma Openings Up to 3-1/2"

## MARLEN UltraPlus™ One-Piece Disposable System For Ileostomies and Colostomies

- With Marlen's Skin Shield™ Adhesive Skin Barrier
- Stoma sizes from 2" (50mm) to 3-1/2" (88mm) and 9 oval sizes from 1" x 1-1/2" (25mm x 38mm) to 2" x 3-1/2" (50mm x 87mm)
- Large pouch size only
- Choice of transparent or opaque outer side



## MARLEN Infant and Pediatric Ileo Pouches

### Infant

- Starter hole is 1/4" (6mm). Cut-To-Fit stoma sizes up to 1-1/4" (32mm)
- With Marlen's Skin Shield™ Adhesive Skin Barrier
- Capacity – 5 oz. (150cc)



### Pediatric

- Starter hole is 1/4" (6mm). Cut-To-Fit stoma sizes up to 1-3/4" (44mm)
- With Marlen's Skin Shield™ Adhesive Skin Barrier
- Capacity – 6-3/4 oz. (200cc)



For more information, contact:



# MARLEN

**Marlen Manufacturing & Development Company**  
 5150 Richmond Road • Bedford, OH 44146 USA  
 Phone: 216.292.7060 • Fax: 216.292.9196  
 E-mail: [info@marlenmfg.com](mailto:info@marlenmfg.com) • [www.marlenmfg.com](http://www.marlenmfg.com)